## **Instrument Assessment Checklist**

Your patients' oral health is your #1 priority. Keeping your instruments "healthy" is just as important. Instruments that are in poor condition compromise patient comfort and safety. Plus, you have to work harder to achieve optimal clinical outcomes.

## **Step 1:** Gather Questionable Instruments:

Consider how long you have had this instrument, whether it's still in use, and if it's performing at an optimal level.

**Step 2:** Check ALL for the Following (any instruments with these or other compromising conditions should be removed from service):

- Bent/broken/misshapen working end
- Corrosion/pitting/rust/discoloration
- Dull cutting edge or point
- Excessive build-up of residue
- Loose joints
- Scratched surface
- Faded markings

Diagnostic Instruments

Diagnostic mstraments.
☐ Scratched mirrors
$\square$ Faded probe markings
$\square$ Dull or bent explorers
$\square$ Bent dressing plier tips
$\square$ Over-compressed dressing pliers
☐ Corrosion
Scalers & Curettes
$\square$ Retipped instruments
$\square$ Dull/chipped cutting edge or point
$\square$ Broken or bent tips
☐ Thin/worn-out blades
☐ Pointed curettes
Ultrasonic Scaling
$\square$ Worn/bent tips
☐ Warped and/or separated stacks
$\square$ Cracked o-ring seals or grips
☐ Water leaking from handle
☐ Bent or distended water delivery tube

Hinged Instrument	ts	
☐ Loose joints		
☐ Corrosion/pitting	, ,	
☐ Rust/discoloration		
•	sshapen working end	
☐ Scratched surface		
	•	
	issing latches g/Rust/discoloration sshapen shell or hinges	
	ized Hu-Friedy representative for more infonts that need to be replaced: (Please list parties)	
Part Code	Description	Qty
		·
		. <u> </u>
		<u> </u>
		. <u>——</u>
		·
		<u> </u>
		· <u></u>
		· —
		·
		<u> </u>
		·
		<del></del>
	-	<u> </u>
		· <u></u>
		· <u></u>
		· <u></u>
		·
		<del></del>
		<u> </u>