For Immediate Release

SATO's PJM RFID Technology Integrated with BIT Consulting's myMEDISET Solution for orthopedic implant vendors based on SAP software

Tokyo, **June 24**, **2015** – SATO (TSE:6287), a global provider of Auto-ID solutions that empower workforces and streamline operations, and BIT Consulting, a leading software solutions provider for the healthcare industry, announced today that SATO's PJM RFID Technology is now available with myMEDISET. BIT Consulting's myMEDISET is a leading solution for medical device loan-kit and consignment management that is based on SAP® software. The myMEDISET 7.0 solution is certified by SAP to run on the SAP NetWeaver® technology platform.

Founded in Freiburg, Germany, BIT Consulting developed myMEDISET for businesses in the medical field. The solution minimizes costs and reduces administrative workloads for managing high-value assets. myMEDISET supports complicated loan-kit operations for orthopedic implant vendors using barcode technology, and with the integration of SATO's PJM RFID Technology, myMEDISET is now an RFID-ready solution.

"Integrating PJM RFID technology with myMEDISET adds real value for our customers," said Michael Knauer, CEO and founder of BIT Consulting. "When a loan kit is dispatched for surgical treatment, it needs to be monitored closely, and employing PJM RFID drastically improves loan processing. The orthopedic industry in Australia has almost universally deployed PJM RFID technology, resulting in vendors and hospitals achieving extremely high levels of efficiency and some of the lowest DPPM rates in the world."

"Many of our customers in Australia, Japan and other regions are using myMEDISET for one reason: it works," said Hiroyuki Konuma, President of SATO Healthcare. "There is nothing more important for our customers when making the jump to RFID. It has to work, and it has to work well. myMEDISET is a proven robust solution used by the world's top

brands, and we are excited to bring our technologies together to benefit customers around the world."

myMEDISET is available with PJM RFID technology via a software upgrade, and is jointly supported by BIT Consulting and SATO. The solution will be exhibited at the Healthcare Information and Management Systems Society show, HIMSS APAC 2015, in Singapore from September 6 to 10, 2015.

About BIT Consulting

BIT Consulting is a world market leader in providing add-on solutions for SAP software for medical device companies, especially in the area of managing complex loan and consignment processes. BIT Consulting's solution myMEDISET has been established over the years as the standard solution for the medical device industry to support efficient management of high-quality assets. Highly customer oriented, BIT Consulting has now been serving customers worldwide for more than 15 years. More information about BIT Consulting can be found at www.bit-consulting.com.

About SATO

SATO (TSE:6287) is a leading global provider of Auto-ID solutions that connect people, goods and information. It serves a diverse range of customers, delivering end-to-end solutions that streamline operations, empower workforces and help customers reduce their environmental impact. For the fiscal year ended March 31, 2015, it reported revenues of JPY 99,831 million (US\$ 910 million*). More information about SATO can be found at www.satoworldwide.com or www.linkedin.com/company/sato-worldwide.

*Conversion is based on an average exchange rate of 1 US Dollar = 109.75 Japanese Yen.

###

SAP, SAP NetWeaver and other SAP products and services mentioned herein as well as their respective logos are trademarks or registered trademarks of SAP SE (or an SAP affiliate company) in Germany and other countries. See http://www.sap.com/corporate-en/legal/copyright/index.epx for additional trademark information and notices. All other product and service names mentioned are the trademarks of their respective companies.

For Media Enquiries:

Daphne Tay
Public Relations Group
SATO Holdings Corporation
daphne.tay@sato-global.com

Tel: +65-9665-2145

SCRUTTINE I